

# ORIVET GENETIC PET CARE MIXED BREED IDENTIFICATION DNA TEST & LIFE PLAN

*Congratulations on making the journey to understanding your dog's ancestry and health risks!*

The Orivet Breed Identification DNA Test and Life Plan helps you discover the secrets buried in your dog's DNA. Understanding your dog's ancestry can assist you and your Veterinarian in creating an effective and personalised training, wellness, nutrition and medical program for your dog.

The greater your knowledge is about your dog's past, the smarter you can be about their future. Dogs - like people - inherit physical and personality traits from breeds in their family tree. They also pass along breed-specific genetic disorders from generation to generation.

The Orivet Breed Identification DNA Test and Life Plan allows you to:

- ✓ Determine the breeds that make up your dog
- ✓ Gain insights into your dog's personality and behaviour
- ✓ Discover genetic health secrets
- ✓ Develop a personalised health and wellness plan
- ✓ Screen hereditary diseases relevant to the breeds that make up your dog
- ✓ Adult weight prediction
- ✓ Improve communication with your veterinarian and other pet care providers



## TERMS AND CONDITIONS:

The Orivet Genetic Pet Care Breed Identification DNA Test provides a genetic insight into the physical and behavioural traits in breeds that have been identified in your dog. The test is designed to be used solely to identify the genetic make up of a mixed-breed dog and no other purpose is intended, authorised, or permitted, including diagnosing diseases or predicting behaviour. Every dog is unique and their temperament and physical traits are the result of many factors, including genetics, environment and training.

Upon receipt of your sample, your dog's DNA is analysed to determine the breeds present in its ancestry with a reasonable degree of certainty. Our procedures are designed to provide reliable and accurate results. Orivet Genetic Pet Care is not responsible for any errors in obtaining the swab sample or for any injuries or loss that may occur as a result.

The purpose of the Breed Identification DNA Test is to identify the genetic makeup of a mixed breed dog. The test is not designed to validate the purity of a purebred dog, and test results should not be relied upon as official certification of your dog's genetic makeup relating to the laws regarding restricted breeds.

The Orivet Genetic Pet Care test is based upon the database of validated breeds. If your dog contains other breeds, it may result in the identification of breeds, or a combination of breeds in the database, related to the breeds in your pet's history.

In the unlikely event that it is not possible to determine breed history or that an error in the analysis occurs, liability by Orivet Genetic Pet Care or related companies and individuals is disclaimed and damages in any event are limited to the payment actually received by Orivet Genetic Pet Care for the specified analysis at issue.

# MIXED BREED IDENTIFICATION DNA TEST & LIFE PLAN

Get to truly know yo♥r best friend



\_\_\_\_\_/\_\_\_\_\_/\_\_\_\_\_  
Date: (Sample Posted)

Pet Name: \_\_\_\_\_

Place Barcode Here

Keep this for your Reference

# IT'S EASY AS 1, 2, 3.

## 1 COLLECT YOUR PET'S DNA

Swab Collection Instructions – view our DNA collection video at [www.orivet.com](http://www.orivet.com)

Please wait at least 30 minutes after a meal or treat (drinking water is fine) to begin your DNA collection. Before collecting the sample, inspect your animal's mouth between cheek and gums for any food material.

**A. COLLECT DNA (EPITHELIAL CELLS)** – Open swab packet carefully and avoid touching the bristle head with your hands at all times. Place the bristle (nylon) head against the inside of the cheek. Pinch the cheek around the bristle head and firmly roll the swab around the pinched area. Pinch and swirl for about 15 seconds. Repeat this step with the remaining swabs.

**B. DRY THE SWABS** – Allow the swabs to air dry for 3-5 minutes. This can be done by standing the swabs up in the allocated holes found in your DNA kit box.

**C. REINSERT DRY SWABS** – After the swabs have dried, place them back into the swab packet. Use the barcode sticker to close the packet, but do not seal it completely. Place a second barcode sticker on the front of the swab packet and keep the remaining barcode for your records (front of this brochure).

## 2 ACTIVATE YOUR KIT ONLINE

Visit [www.orivet.com/activate/kit](http://www.orivet.com/activate/kit) and follow the prompts.



## 3 MAIL YOUR SAMPLE

Please post samples back to us using the self-addressed postage-paid envelop provided. Samples should be stable at room temperature

**YOU WILL RECEIVE AN EMAIL TO NOTIFY YOU WHEN YOUR SAMPLE ARRIVES**



1



2



3



**Orivet**  
Genetic Pet Care